

# SBFSEM-tools

Sara Patterson

Neitz Lab, University of Washington

December 17, 2017

Goal: basic analysis support requiring minimal setup and background in computers, math, etc.

Current capabilities:

- ▶ Import and parse data into Matlab through Viking's OData service
- ▶ Parse data into high-level data structures (Neuron, etc)
- ▶ Count up synapses and condense synapses spanning multiple slices
- ▶ Summarize synapses with statistics and plots
- ▶ Arbitrary geometry rendering
- ▶ Basic network analysis

# Install

First download or clone [SBFSEM-tools](#).

Make sure SBFSEM-tools is added to your MATLAB [path](#) like so:

```
1 addpath(genpath('C:\...\sbfsem-tools'));
```

If you already have JSONLab installed, make sure

```
1 which loadjson
```

returns the version in sbfsem-tools. Otherwise, you might get some errors.

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

Appendix

Links

References

Tulip Import

NeuronApp

# Neuron

The Neuron class is the basic representation of a cell exported from Viking. To create a Neuron:

```
1 % cellName = Neuron(cellID, 'source');  
2 c207 = Neuron(207, 'temporal');
```

Note: Source can be 'inferior', 'temporal' or 'rc1'. Also, abbreviating to 'i', 't' and 'r' works as well.

Open the Neuron class in the user interface

```
1 NeuronApp(c207);
```

To update a neuron in the workspace:

```
1 c207.update();
```

# Neuron

## Properties

Neuron has the following publicly accessible properties:

- ▶ data
- ▶ viking - Structure of the neuron's info from viking
- ▶ nodes - Table of all annotations
- ▶ edges - Table of links between annotations
- ▶ volumeScale - Units are nm/pix for X,Y and nm/section for Z
- ▶ synapses - Table of each child structure.
- ▶ geometries - Table of closed curve geometries
- ▶ analysis - Containers.map of analyses
- ▶ lastModified - Last update of neuron from OData

[SBFSEM-tools](#)

[Sara Patterson](#)

[SBFSEM-tools](#)

[Install](#)

[Neuron Class](#)

[Analysis](#)

[Renders](#)

[Export](#)

[Appendix](#)

[Links](#)

[References](#)

[Tulip Import](#)

[NeuronApp](#)

# NeuronGroup

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

Appendix

Links

References

Tulip Import


NeuronApp

Create a single data object to hold related Neurons. Note inputs can be ID numbers or existing Neurons.

```
1 h1hc = sbfsem.NeuronGroup([28, 447, 619]);  
2 % Add a neuron to an existing group  
3 h1hc.add(4568);  
4 % Remove a neuron  
5 h1hc.remove(4568);
```

To plot the somas of all Neurons in the NeuronGroup:

```
1 h1hc.somaPlot();  
2 h1hc.somaPlot('addLabel',true); % Label with ID  
3 h1hc.somaPlot('ax',gca); % Add to existing axis  
4 % Two methods for controlling plot color:  
5 h1hc.somaPlot('Color',[0 0.8 0.3]);  
6 h1hc.setPlotColor([0 0.8 0.3]); h1hc.somaPlot;
```



Mosaic of H1 and H2 somas

The NeuronAnalysis class helps keep population data organized by managing input parameters and results of common analyses. To create a new analysis, subclass NeuronAnalysis and edit the doAnalysis and visualize methods.

See Tutorial.m for information on these existing classes:

- ▶ **DendriticFieldHull** - uses convex hull to estimate dendritic field area, includes methods for removing axons prior to analysis.
- ▶ **PrimaryDendriteDiameter** - returns the median dendrite diameter at a given distance from the soma.



```
1 c6800 = sbfsem.Neuron(6800, 'i');  
2 sbfsem.ui.StratificationView(c6800);  
3 c207 = sbfsem.Neuron(207, 'i');  
4 sbfsem.ui.SomaDistanceView(c207);  
5 c4781 = sbfsem.Neuron(4781, 'i');  
6 sbfsem.ui.NodeView(c4781);  
7 sbfsem.ui.NodeView(c6800);
```



# IPL Boundary Surfaces

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

Appendix


Links

References

Tulip Import

NeuronApp

```
1 inl = sbfsem.core.INLBoundary('i');  
2 % Update marker locations from OData  
3 inl.refresh();  
4 % Create the surface  
5 inl.doAnalysis();  
6 % Plot the surface  
7 plot(inl);
```



# XY Offset

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

Appendix

Links


References

Tulip Import

NeuronApp

The function `xyRegistration.m` calculates the XY offset through a range of Z sections and outputs statistics on the offsets relative to the most sclerad section input.

```
1 % S xyRegistration(source, [minZ maxZ], plotFlag);  
2 S = xyRegistration('i', [1284 1309], true);
```



Next on the todo list is a systematic way of applying transforms

VikingPlot generates 3D models by fitting the data to three stereotyped geometries - sphere, cylinder and cone. I'm focusing my efforts on the opposite approach - rendering the structures 'as is'. By not fitting the data, the renders will be far more accurate. The tradeoff is that this accuracy applies both to the neuron's morphology and the small discrepancies in image registration.

# Closed Curve Renders

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

**Renders**

Export

Appendix

Links



References

Tulip Import

NeuronApp

Closed-curve structures can be rendered into volumes:

```
1 c2542 = sbfsem.Neuron(2542, 'i');  
2 Imcone = sbfsem.render.ClosedCurve(c2542);
```




See the Tutorial.m file for more information on render colors, lighting and materials.







# VikingPlot Comparison



Efficient: This render took 32.77 sec on my laptop.

[SBFSEM-tools](#)

[Sara Patterson](#)

[SBFSEM-tools](#)

[Install](#)

[Neuron Class](#)

[Analysis](#)

**[Renders](#)**

[Export](#)

[Appendix](#)

[Links](#)

[References](#)

[Tulip Import](#)


[NeuronApp](#)

# Disc Renders

## Stacked Cylinders method

The Disc Render uses a similar method to the Closed Curve renders. While the result does have 3D lighting, it needs further improvement before being a viable alternative to VikingPlot or rotated cylinders.

```
1 c1893 = sbfsem.Neuron(1893, 'i');  
2 rodBC = sbfsem.render.Disc(1893);
```



[SBFSEM-tools](#)

[Sara Patterson](#)

[SBFSEM-tools](#)

[Install](#)

[Neuron Class](#)

[Analysis](#)

[Renders](#)

[Export](#)

[Appendix](#)

[Links](#)

[References](#)

[Tulip Import](#)


[NeuronApp](#)

# ImageStackApp

The ImageStackApp was built to scan through a series of images exported using Viking's Export Frames option (although this will work with any folder of images). The app will import all images in a given folder. If the filenames contain numbers (Export Frames appends a frame number automatically), these will be used to order the images. Use the buttons (or the left/right arrow keys) to browse through the images.

```
1 ImageStackApp( 'C:\...\foldername' );
```

To crop the images, use Process Image  
→ Crop and draw a rectangle on the image. The change will not be applied until moving to another image.



SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

Appendix

Links

References

Tulip Import

NeuronApp

# Appendix

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

**Appendix**

Links

References

Tulip Import

NeuronApp

Here's some information and methods that are less essential:

1. Resources
2. Links
3. Tulip import method
4. Old version of NeuronApp

The SBFSEM-tools [repository](#) can be found on Github.

The software used for annotations is [Viking](#), developed by Jamie Anderson and the Marc Lab at University of Utah. Useful free programs involved in these analyses:

- ▶ [Tulip](#) supports graph visualization. The documentation for Tulip's Python package can be found [here](#).
- ▶ [Blender](#) for 3D renders of neurons.

SBFSEM-tools was developed in the [Neitz lab](#) at University of Washington.

## 3D rendering

- ▶ Lorensen & Cline (1987) *Computer Graphics*, 21(4)
- ▶ 3D Math Primer For Graphics and Game Development by Fletcher Dunn and Ian Parberry
- ▶ Blender Cheatsheet

## OData

- ▶ Microsoft Developer OData tutorials

## Data Structures

- ▶ BaseCS
- ▶ Data structure visualization

# Import

## Step One: Tulip

Open a cell in Tulip and then open the Python command line. It's on the bottom toolbar.

Set the file name\* and file path:

```
1 outputFile = "C:\...\c207.json"
```

Then run these two lines:

```
1 params = tlp.getDefaultPluginParameters( 'JSON  
Export ', graph)  
2 success = tlp.exportGraph( 'JSON Export ', graph  
, outputFile, params)
```

# Import

## Step One: Tulip

Tulip 4.19 [Tulip] - C:\Users\sarag\Google Drive\NEITZ Lab\Electron Microscopy Rn\SBFSEM-tools\demo\sb800.bpx  
File Edit Window Help

**Elements**

Show **Items** 2 filtered by **no selection**

matching given a pattern **id** in **layer 2** column

show columns matching given a pattern

localid	localidNew	switchedActivation	CER
0	27188	0.00017500000007	2.0
1	27187	0.00010000000000	2.0
2	27188	0.00010000000007	2.0
3	27188	0.00010000000007	2.0
4	27188	0.00010000000007	2.0
5	27187	0.00010000000000	2.0
6	27187	0.00010000000000	2.0
7	27187	0.00010000000000	2.0
8	27187	0.00010000000000	2.0
9	27187	0.00010000000000	2.0
10	27187	0.00010000000000	2.0
11	27187	0.00010000000000	2.0
12	27187	0.00010000000000	2.0
13	27187	0.00010000000000	2.0
14	27188	0.00010000000007	2.0
15	27188	0.00010000000007	2.0
16	27188	0.00010000000007	2.0

**Selected graph:** TLP import - Name: C:\Users\sarag\Google Drive\NEITZ Lab\Electron Microscopy Rn\SBFSEM-tools\demo\sb800.bpx

```
Python 3.7.3 64-bit, Apr 10 2017, 23:31:20 [AMD v.1500 32 bit (Intel)] on win32
# Use Ctrl+L to show dynamic info-completion dialog
>>> outfile = "C:\Users\sarag\Google Drive\NEITZ Lab\Electron Microscopy Rn\SBFSEM-tools\demo\sb800.bpx"
>>> param = "TLP import - Name: C:\Users\sarag\Google Drive\NEITZ Lab\Electron Microscopy Rn\SBFSEM-tools\demo\sb800.bpx"
>>> success = tlp.exportGraph('SBFSEM Export', graph, outfile, param)
>>>
```

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

Appendix

Links

References

Tulip Import

NeuronApp



# Alternative Streamlined Import

## Step One - Tulip

If you're comfortable with Python, you can skip opening Tulip's UI entirely.

First get the Tulip modules:

```
1 $ pip install tulip-python
```

This will allow you to export .t1p (or compressed .t1p.gz) files to JSON from a command line without Tulip's UI.

```
1 from tulip import tlp
2 graph = tlp.loadGraph("C:\...\morph-207.t1p")
3 outputFile = "C:\...\c207.json"
4 params = tlp.getDefaultPluginParameters('JSON
   Export', graph)
5 success = tlp.exportGraph('JSON Export', graph,
   outputFile, params)
```

# Import

## Step Two - MATLAB

Load in the JSON file and create a Neuron object:

```
1 c207 = Neuron( 'c207.json' );
```

A dialog box will ask for the cell number and source (temporal, inferior, rc1). To avoid that, include them while creating the Neuron object:

```
1 % output = Neuron(filename, cellNumber, source);  
2 c207 = Neuron( 'c207.json', 207, 'temporal' );
```

To update the underlying data for an existing Neuron:

```
1 c207.updateData( 'c207.json' );
```

Open up the user interface:

```
1 c207.openUI;
```

# Cell Info Panel

This entire panel is designed with a future directory class in mind. As of now, you don't need to save each Neuron, so only set these if you have a specific reason for doing so.

- ▶ If known, the cell type and subtype will be helpful for connectivity analysis.
- ▶ The other properties will eventually be used for cell queries but don't have much use yet.

After changing any of the attributes on the Cell Info Panel, make sure to press the [Add to cell data] button. This will make sure your changes are reflected next time you open the UI.

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export

Appendix

Links

References


Tulip Import

**NeuronApp**


# 3d plot

## Components

You can add and remove each synapse type, the soma node and the skeleton independently using the checkboxes. Rotate the plot with the elevation and azimuth sliders.



All amacrine cell synapses



Semilunar synapses & skeleton


# Histograms

There are two histograms.

These plot synapse count as a function of:

- ▶ Distance from soma
- ▶ Section number (z-axis)

Use the Synapse Table to edit the number of bins. Add the cell skeleton to see dendrite stratification.




Z-axis synapse distribution for a putative AII AC. This synapse asymmetry isn't news but is still good to see.

# Connectivity

To add connectivity data, save a network map in Tulip as a JSON file, as described above. Then in Matlab:

```
1 c207.addConnectivity('c207hops.json');
```



The connectivity matrix is weighted by the number of unique synapses between two cells (dark blue for 0 synapses). Directed synapses will only register a contact from the pre → post-synaptic neuron.

Connectivity for the All AC.

The network data is split into edges and nodes, each with their own table. The tables can be exported to Excel as .csv or as a text file from the UI menu bar.

To print a easily readable table to the command line use the `networkTable` function

```
1 networkTable(c207);
```

# Blender Renders and Cell Fills

SBFSEM-tools

Sara Patterson

SBFSEM-tools

Install

Neuron Class

Analysis

Renders

Export



Appendix

Links

References

Tulip Import

**NeuronApp**



For now, set the `renderDir` in `getFilepaths.m` to the file your images are saved into. The UI find the images if their filename includes the letter 'c' followed by cell number (like 'c207'). I hope to improve this at some point.

This isn't limited to renders and could include whatever images and diagrams are helpful. For example, I can compare my ON-smooth cell reconstructions and cell fills.